

# EasyOne® Air

*Simple. Certain. Proven.*



## Innovators for Optimal Lung Care

EasyOne Air is the new, all-in-one revolution in lung function testing. Built for healthcare providers large and small, EasyOne Air puts the power of the most advanced spirometry testing tools in the palm of your hand.

Revolutionary TrueFlow ultrasound technology means unprecedented accuracy for reliable diagnosis and trusted treatment. Equipped with ready-to-use EMR connectivity and a user-friendly interface, the EasyOne Air makes the best in lung care available anywhere.

Learn more about EasyOne Air at [nddmed.com](http://nddmed.com)

300 Brickstone Square, Suite 604 Andover, MA 01810 | (978) 470-0923 | [customerservice@nddmed.com](mailto:customerservice@nddmed.com)



# Simple. Certain. Proven.

Rechargeable battery

Wireless connectivity to the pc via Bluetooth

Easy navigation

- Quickly assesses test quality with automated QC grades and messaging

Large Touch Screen

Incentive screens available with EasyOne Connect

- Instant interpretation
- EasyOne Connect for seamless EMR integration
- Internal storage for up to 10,000 tests with SD card
- 3 year exchange warranty
- Automatic calibration
- Immediate test quality feedback

**TrueFlow**  
*makes the difference*

The original ultrasonic flow measurement is highly accurate in all flow ranges, independent of gas composition, pressure, temperature and humidity. ndd TrueFlow is a hygienic and resistance-free solution that does not require calibration during its lifetime.

**EasyOne Connect**  
*intelligent workflows*

EasyOne Connect software offers one database and one platform for all EasyOne point-of-care-solutions. Data management has never been easier.

## Specifications

**Test Types** FVC, FVL, SVC, MVV and Pre/Post Testing

**Languages** English, French, German, Spanish

**Printing Options** Direct to printer or with EasyOne Connect software

**EMR connectivity** EasyOne Connect

**Dimensions** 3.4x6.1x1.4", 13 OZ

**Meets ATS/ERS 2005, OSHA and SSD requirements**

## Part Numbers

**2500-2A** EasyOne Air spirometry system

**5050-50** EasyOne Flow Tube box of 50

**5050-200** EasyOne Flow Tube box of 200