## OPTICLUX 10-2 C T1 (Wood-Light) D16071000



ENGINEER OF LIGHT.











fitted with LED - Energy efficiency category A+

color rendering index (CRI) > 95 colour temperatur 6500 K

work equipment electronic converter mains lead direct power supply

lens; dimensions Ø 160 mm lens; diopters 3,5

luminaire body polyamide (PA) weight (net) 9.1 kg (20 lb)

Power consumption 12 W

usage keypad at luminaire head

balance of articulation joints spring class of protection I

Accessory adapter for suspended ceiling norms EN 60601-1, EN 60601-2-41

special features Rotary feedthrough (360° Turnover axis),

scratch-proof plastic lens, continuously dimmal

switchable Wood-Light

Fastening ceiling mount, ceiling tube can be cut to length

(40-116cm)

## light distribution





