

# The FastPack® IP System Specifications

## Analyzer Specs

- Technology: Paramagnetic Particle Chemiluminescence Immunoassay
- Size: 13 x 9 x 12 inches (33 x 23 x 30 cm)
- Weight: 28.64 lbs (12.99 kg)
- Electrical Requirements: 100-240  $\pm$  10% VAC, 1.9A
- Power Consumption: 100 Watts (maximum)
- Ambient Operating Temperature: 15°C (59°F) to 32°C (90°F)
- Operating Humidity: 10% to 80% relative humidity
- Operating Altitude:  $\leq$  2000 m (6562 ft.)
- Display Output: 2-line alphanumerical display



## Printer Specs

- Size: 7.1 x 5.1 x 7.3 inches (18 x 13 x 18.5 cm)
- Weight: 6 lbs (2.7 kg)
- Electrical Requirements: 90-264 V AC, 50/60 Hz (input); 24V DC  $\pm$  5%, 2.0A (output)

The **FastPack® IP System** — Diagnostic technology that gives you the power to provide reliable results in real-time.

**QUALIGEN DIAGNOSTICS**  
2042 Corte Del Nogal, Suite B  
Carlsbad, California 92011  
877.709.2169 option 3  
FastPack@qualigeninc.com

