TRIANGO 80-1 C D16138000



ENGINEER OF LIGHT.











Illuminancefield @ 100 cm Ø 18 cm (7.1 in) illuminance @ 100 cm 80 000 lx color rendering index (CRI) > 95 4500 K color temperature depth of illumination L1+L2 1.24 m (48.8 in) work equipment electronic converter mains lead direct power supply

luminaire body polycarbonate & acrylonitrile butadiene styrene

copolymer (PC & ABS)

weight (net) 14.0 kg (30.9 lbs)

30 W power consumption

usage rotary switch on the luminaire head

single-pole secondary, continuously dimmable

balance of articulation joints spring (weight compensation)

class of protection

accessory sterilizable handle, adapter for suspended ceiling AAMI/IEC 60601-1; CAN/CSA-C22.2 60601-1 norms IEC 60601-2-41; CAN/CSA-C22.2 60601-2-41

AMD 1/IEC 60601-1-6; CAN/CSA-C22.2 60601-1-6

special features sterilizable grip, continuously dimmable, maintenance

fastening ceiling mount, when mounting ceiling tube could be

cut to length 25-120 cm (9.8-47.2 in)

system of protection **IP 43**

order number D16138000 - pure white

light distribution





