Examination luminaire

VISIANO 10-1 P S10 D15757150



rendering index (CRI) > 93

color temperature work 4400 K

equipment electronic converter connected load 100-240 V; 50/60 Hz

mains lead 1.8 m (5.9 ft) hospital grade plug

luminaire body PMMA weight (net) 1.1 kg (2.4 lbs)

power consumption 13 W

usage rocker switch (I/0) at the arm

class of protection II

accessory rail clamp

norms EN 60601-1, EN 60601-2-41

special features positioning handle certifications ETL/cETL system of protection IP 20

order number D15757150 - pure white

Derungs D

light distribution



