## **Examination luminaire**

## VISIANO 20-2 P TX D15461150



ENGINEER OF LIGHT.



fitted with 1 x LED module
luminous power 60 000 Lux / 0.5 m
Illuminancefield 0.5 m Ø 21 cm (8.3")

Color Rendering Index > 95

(CRI)

colour temperatur 3500 K & 4500 K
work equipment electronic converter
connected load 100-240 V; 50/60 Hz

mains lead 1.8 m (5.9 ft) Hospital Grade Plug

luminaire body polycarbonate & acrylonitrile butadiene styrene copolymer

(PC & ABS)

**weight (net)** 2.0 kg (4.4 lbs)

power consumption 21 W

**usage** keypad at luminaire head

balance of articulation spring

joints

. .

class of protection

accessory included rail mount bracket with removal protection

**norms** EN 60601-1, EN 60601-2-41

special features 2 color shades, 4 dimming levels, antimicrobal, Rotary

feedthrough (360° Turnover axis)

 $\begin{array}{ll} \text{safety} & \text{ETL/cETL} \\ \text{system of protection} & \text{IP 20} \end{array}$ 

order number D15461150 - pure white

## light distribution

Derungs D





