



fitted with	1 x LED module
luminous power	60 000 Lux / 0.5 m
Illuminancefield 0.5 m	Ø 21 cm (8.3")
Color Rendering Index (CRI)	> 95
colour temperatur	3500 K & 4500 K
work equipment	electronic converter
connected load	100-240 V; 50/60 Hz
mains lead	1.8 m (5.9 ft) Hospital Grade Plug
luminaire body	polycarbonate & acrylonitrile butadiene styrene copolymer (PC & ABS)
weight (net)	2.0 kg (4.4 lbs)
power consumption	21 W
usage	keypad at luminaire head
balance of articulation	spring
joints	
class of protection	II
accessory	included wall bracket
norms	EN 60601-1, EN 60601-2-41
special features	2 color shades, 4 dimming levels, antimicrobial, Rotary feedthrough (360° Turnover axis)
safety	ETL/cETL
system of protection	IP 20
order number	D15461120 - pure white



light distribution

